



**Malvern
Panalytical**

Biophysical characterization of SMALPs and nanodiscs


Verna Frasca, Ph.D.
Malvern Panalytical
Verna.Frasca@Malvern.com

Malvern Panalytical



Provider of industry leading analytical instrumentation for the characterization of

Molecules



Particles



Materials



Along with the expertise and understanding of how

Molecular
Properties

Control


Particle Size,
Shape, &
Interactions


Drives


Product
Performance


Malvern Panalytical Biosciences Group Solutions and Instrumentation for SMALP and nanodisc characterization



| Zetasizer (0.6 nm – 10 μ m) ★ | |
|---|---|
|  | Light Scattering |
| | Stability Screening & Aggregate Detection |
| | PSD, R_H , T_{Agg} , k_D , B_{22} , & Z_{Eff} |

| OMNISEC/MALS ★ | |
|---|--|
|  | Advanced Detection SEC |
| | Stability Screening |
| | M_W , R_H , %Purity, & HMW Detection |

| MicroCal PEAQ DSC ★ | |
|---|-----------------------------------|
|  | Differential Scanning Calorimetry |
| | Stability Screening |
| | T_M & ΔH |

| MicroCal PEAQ ITC ★ | |
|---|---|
|  | Isothermal Calorimetry |
| | Activity Screening |
| | K_D , ΔH , ΔS , & n |

- **Driving advances in biophysical characterization** by engaging and collaborating with researchers to improve understanding of biomolecular interactions and to accelerate development of more effective drugs
- From molecular interactions to successful drug products

Zetasizer Nano & APS


Submicron Light Scattering

- **Size, Charge, and Interaction Parameters**
- Proteins, peptides, biopolymers, and nanoparticles
- Screen for aggregation propensity
- Zeta potential
- Rapid aggregate assessment
- **Screen For Colloidal Stability**
- Multivariate data sets for stability and aggregation metrics from ultra-low volume/concentration assays


SLS | DLS | ELS




R_H | R_S | k_D
 Z_{Eff} | B_{22} | T_{Agg}
Distribution



Reversible Self-Association



Electrostatic Repulsion



Stability Profile

IgG Stability Profiles For Buffer 1 And Buffer 2

| Property | Buffer 1 | Buffer 2 |
|---|----------|----------|
| k_D (ml/g) | -5.2 | 31.9 |
| B_{22} ($\times 10^5$ ml mole/g ²) | -1.5 | 127.5 |
| Z_{Eff} | 0.7 | 4.3 |
| T_m (C) | 56 | > 56 |
| T_{Agg} (C) | 66 | > 66 |
| R_S (nm) | 5.8 | 5.7 |
| SubQ Limit (mg/ml) | 137 | 148 |

Probing molecular interactions of poly(styrene-co-maleic acid) with lipid matrix models to interpret the therapeutic potential of the co-polymer



Particle size and zeta potential measurements with Zetasizer NanoZS DLS

| System ^a | Hydrodynamic diameter ^b (nm) | Zeta potential ^b (mV) |
|---------------------|---|----------------------------------|
| A | 1453 ± 21.5 | -3.62 ± 0.42 |
| B | 1746 ± 31.8 | -8.54 ± 0.72 |
| C | 1501 ± 24.2 | -4.18 ± 0.47 |
| D | 1739 ± 28.7 | -8.26 ± 0.57 |
| E | 1420 ± 38.1 | +12.9 ± 1.12 |
| F | 1579 ± 25.4 | -8.72 ± 0.50 |
| G | 1507 ± 42.9 | -17.1 ± 1.54 |
| H | 1852 ± 36.3 | -26.5 ± 1.06 |
| | 194.8 ± 1.4 | |

^a Compositions of the systems are as follows: A, DSPC; B, DSPC + SMA; C, DSPC + CHOL; D, DSPC + CHOL + SMA; E, DSPC + DODAB + CHOL; F, DSPC + DODAB + CHOL + SMA; G, DSPC + DCP + CHOL; H, DSPC + DCP + CHOL + SMA.

^b Mean ± SD (n = 3).

- Samples: Neat and SMA-incorporating MLVs at pH 7.0
- The observed slight increase in the mean hydrodynamic diameter of the MLVs in presence of SMA was due to incorporation of the co-polymer within the bilayers
- Zeta potential measurements revealed significant change in the vesicle surface charge in presence of the anionic co-polymer.

Banerjee et. al, Biochim Biophys Acta, 1818, 537-550 (2012)

OMNISEC/MALS

Advanced Detection & MALS SEC

- **Absolute M_w , oligomeric distribution, %Purity, and size from a single injection**
- Multi detection SEC includes UV, RI, MALS, and DLS, as well as DSV for intrinsic viscosity
- Increased light scattering sensitivity for detecting trace amounts of HMW species
- **Define Oligomeric Baseline**
- Measure %Purity & distribution
- Track aggregation & HMW species generation

SEC/LS | MALS | DLS | DSV




R_H | M_w | R_g
 C | $[\eta]$ | Shape
 %Purity | Distribution
 HMW Species Detection


Denosumab (Prolia® and Xgeva®): Innovator and Biosimilar

Stressed – Incubation at 30°C

| Innovator | Biosimilar |
|-------------|-----------------|
| 99% monomer | 97% monomer |
| 1% dimer | 1.5% dimer |
| | 1.5% aggregates |



SMALP analysis (SMA distribution) – OMNISEC data




UV area: 192 mVmL

RI, RALS, DP, A_{260}

UV area [1]: 99 mVmL

→ Nearly 50% of SMA is free in solution


UV area [2]: 93 mVmL



Johannes Klingler, Malvern Panalytical webinar, 2017



SMALP analysis (SMA distribution) – OMNISEC data



C_{SMA} , C_{DMPC} , C_{total}

→ The second peak is virtually lipid-free




Johannes Klingler, Malvern Panalytical webinar, 2017



© Malvern Panalytical 2018

SMALP analysis (SMA distribution)



C_{SMA}, C_{DMPC}

→ SMA is more abundant in the right flank of the SMALP peak



Johannes Klingler, Malvern Panalytical webinar, 2017



© Malvern Panalytical 2018

SMALP analysis (SMALP size and composition)

– OMNISEC data



M_w : 161 kg mol⁻¹

M_w/M_n : 1.05

$[\eta]$: 0.049

R_h : 5.0 nm

$X_{SMA/DMPC}$: 0.40

C_{SMA} , C_{DMPC} , C_{total} , R_h , $[\eta]$,
 M_w

→ SMALP masses and sizes are narrowly distributed
→ SMA contributes 40% to the total mass of SMALPs

Summary and conclusions: OMNISEC



- Nearly 50% of SMA is free in solution
- SMALPs show relatively narrow size and mass distributions
- SMA contributes a significant amount to the SMALP mass
- Combining data on hydrodynamic size and mass confirms disc shape of SMALPs


MicroCal PEAQ DSC

Differential Scanning Calorimetry

- **Conformational Stability**
- Sensitive to biopolymer domain transitions
- Phase transition reversibility
- T_m correlated with aggregation propensity, formulation stability, and shelf life
- **Screen For Structural Stability**
- High quality thermal phase transition data
- Gold standard, label-free, universal tool for studying thermal stability



Conformational stability




Multiple descriptors from thermal stability profile for fingerprinting



T_m | T_{Onset} | $\Delta H_{U...}$

HOS Similarity




Data for Regulatory filing

Probing molecular interactions of poly(styrene-co-maleic acid) with lipid matrix models to interpret the therapeutic potential of the co-polymer



Thermal stability using MicroCal VP-DSC





- Study the thermotropic phase behavior of liposomes from which the molecular interactions between the co-polymer and phospholipids were quantitatively probed
- Representative DSC heating scans of
 - A, DSPC; B, DSPC+SMA;
 - C, DSPC+CHOL; D, DSPC+CHOL+SMA;
 - E, DSPC+DODAB+CHOL; F, DSPC+DODAB+CHOL+SMA;
 - G, DSPC+DCP+CHOL; H, DSPC+DCP+CHOL+SMA
- Multilamellar vesicles prepared in hepes buffered saline (10mM HEPES+150 mM NaCl, pH 7.0).
- Changes in T_M , enthalpy and peak broadening provide insights into the mechanism and interaction between SMA and the MLVs

MicroCal PEAQ ITC




Isothermal Titration Calorimetry

- **Affinity and Function**
- Measures any interactions between two biomolecules in solution or suspension, including nanoparticles and nanodiscs
- Micellization/demicellization
- **Confirm Function**
- High quality affinity data and stoichiometry
- Gold standard, label-free, universal tool for studying biomolecular interactions



FhuA interactions in a detergent-free nanodisc environment

Binding affinity, stoichiometry and thermodynamics using MicroCal VP-ITC



- Binding affinity between Nd-FhuA and TonB32–239. The ITC thermograms show the interaction of Nd-FhuA with TonB in the presence (A) or absence of ferricrocin (B).

Binding affinities and thermodynamics between Nd-FhuA, TonB and ColM.

| Titrant | Cell | N | K _D (nM) | ΔH (cal/mol) | ΔS (cal/mol/deg) |
|-----------|--------------|----------------|---------------------|---------------|------------------|
| TonB | Nd-FhuA apo | - | - | - | - |
| TonB | Nd-FhuA-Fc | 0.977 ± 0.0100 | 200.4 ± 29.3 | -9186 ± 133.4 | -0.164 |
| Colicin M | Nd-FhuA apo | 0.906 ± 0.0016 | 3.48 ± 1.09 | 7969 ± 33.3 | 65.4 |
| Colicin M | Nd-FhuA-Fc | - | - | - | - |
| TonB | Nd-FhuA-ColM | - | - | - | - |

The estimated errors are based on a χ^2 minimized fit of the experimental data to a single-site binding model using Origin 7.0 software (OriginLab).

Summary




- ITC, DSC, DLS, and Multi-detection SEC contribute important information about SMALPs and nanodisc structure, mechanism, and interactions
- Use in conjunction with other techniques for complete biophysical characterization and development of new methods and products
 - SPR
 - FTIR
 - NMR
 - Mass spec
 - CD
 - Analytical ultracentrifugation
 - Microscopy

Malvern Panalytical – Solutions and Instrumentation for Bioscience




MicroCal PEAQ ITC ★




| |
|---|
| Isothermal Calorimetry |
| Activity Screening |
| K_D , ΔH , ΔS , & n |

Viscosizer TD (0.9 – 120 cP)




| |
|--|
| Taylor Dispersion Analysis & Viscosity |
| Stability Screening |
| R_S , k_D , & η |

Morphology G3-ID (1 – 1000 μm)




| |
|--------------------------------------|
| Digital Imaging & Raman Spectroscopy |
| Aggregate Quantification & ID |
| Size, Morphology, & Chemical ID |

OMNISEC/MALS ★




| |
|--|
| Advanced Detection SEC |
| Stability Screening |
| M_W , R_H , %Purity, & HMW Detection |

Empyrean




| |
|----------------------------------|
| X-Ray Diffraction |
| Polymorph Assessment & Screening |
| Excipient & API Form |

Mastersizer 3000 (0.01 – 3500 μm)




| |
|----------------------------|
| Laser Diffraction |
| Quality Control |
| Particle Size Distribution |

MicroCal PEAQ DSC ★




| |
|-----------------------------------|
| Differential Scanning Calorimetry |
| Stability Screening |
| T_M & ΔH |

NanoSight (40 nm – 1 μm)




| |
|--|
| Nanoparticle Tracking Analysis |
| Aggregate Quantification |
| PSD, R_H & $S_{\mu P}$ Concentration |

Kinexus




| |
|---------------------------|
| Rotational Rheometry |
| Quality Control |
| η , G' , and G'' |

Zetasizer (0.6 nm – 10 μm) ★




| |
|---|
| Light Scattering |
| Stability Screening & Aggregate Detection |
| PSD, R_H , T_{Agg} , k_D , B_{22} , & Z_{Eff} |

Archimedes (250 nm – 5 μm)



| |
|--|
| Resonant Mass Measurement |
| Aggregate Quantification & ID |
| PSD, R_m , m_B , & $S_{\mu P}$ Concentration |

Zetium



| |
|-------------------------------------|
| X-Ray Fluorescence |
| Quality Control |
| Residual Catalyst & Metal Detection |

★ SMALP and nanodisc characterization